

# X-Bolt Pro-X1™ Trochanteric Nail

**X-Bolt**  
Trauma Orthopedics

**VALUE+**  
ANALYSIS BRIEF



# VALUE+

ANALYSIS BRIEF



# CONTENTS

BACKGROUND: HIP FRACTURES	2
EPIDEMIOLOGY	3
THE RATIONALE FOR OPERATIVE HIP FRACTURE CARE	3
ECONOMIC BURDEN	4
ECONOMIC IMPACT OF REOPERATION	4
MAJOR CAUSES OF REOPERATION	4
UNMET NEED / SOLUTION	5
FEMORAL HEAD CUT-OUT	6
DISTAL STRESS MODULATION	7
FASTER OPERATIVE TIME	7
ECONOMIC VALUE: THE X-BOLT® PRO-X1™ SOLUTION	8
REFERENCES	9

# BACKGROUND: HIP FRACTURES

Hip fractures present a significant global challenge. As the population ages, osteoporotic fractures become increasingly common, and are escalating costs to health care services worldwide.

1-in-6 white women will suffer a hip fracture in their lifetime. The 1-year mortality rate after a hip fracture is as high as 30%. A large proportion of patients who survive beyond this often require prolonged supportive care.<sup>1,2</sup>


In the UK, nearly 2% of all acute hospital beds are occupied by hip fracture patients and this one injury carries a total cost equivalent of approximately 1% of the whole NHS budget.<sup>3</sup>

## Proximal Femoral (Hip) Fractures (AO/OTA classifications)




# EPIDEMIOLOGY

80% of hip fractures occur in those over the age of 75 years (median age 84 years).<sup>3-5</sup> Approximately 300,000 hip fractures occur each year in the United States.<sup>6</sup> Over 1 million hip fractures occur between the EU27 and North America each year.<sup>6,7</sup> Extracapsular fractures account for approximately 50% of all hip fractures recorded and are generally 'fixed' with a nail or plate. Intracapsular fractures are generally 'replaced' with a hemiarthroplasty or total hip replacement.<sup>1-5</sup>



## THE RATIONALE FOR OPERATIVE HIP FRACTURE CARE

The key driver for operative management is the prevention of pre-terminal decline. Failure to recover some form of mobility renders these already highly vulnerable patients, more liable to develop venous thromboembolism, chest infections and pressure sores. Problems with personal hygiene, independence and activities of daily living rise exponentially with an unfixed fracture. Pain is also a key indication for surgery; from a humanitarian perspective, operative fixation is akin to a palliative procedure to allow comfort for patients in their last days. Hip surgery should be performed urgently, noting that pain and immobility are key indications for surgery. The goal of surgery is to restore mobility and to ameliorate any decline in the patient's independence and ability to function.



# ECONOMIC BURDEN

An analysis of the 2014 Medicare database<sup>8</sup> showed the estimated cost per patient within acute care intertrochanteric hip fractures is **\$52,512** in the first 365 days (\$44,135 in first 90 days). Another paper with data from Medicare 2013-2015 shows one year costs for an intertrochanteric fracture without a reoperation was **\$55,701**.


In Europe<sup>7</sup>, the one-year costs for hip fractures are also substantial, averaging **\$41,832** per patient episode.

EUROPE	Hip Fracture numbers	Annual cost per episode
France	73,630	\$39,717
Germany	129,849	\$36,995
Italy	90,539	\$49,642
Netherlands	13,395	\$47,119
Spain	40,473	\$44,409
Sweden	20,280	\$45,870
UK	79,243	\$29,072
<b>AVERAGE</b>		<b>\$41,832</b>

# ECONOMIC IMPACT OF REOPERATION

Data from Medicare between 2013 and 2015 shows that the additional one year cost for those needing a reoperation following intertrochanteric fracture fixation was **\$38,029** versus those patients who did not have a reoperation (\$93,730 vs \$55,701).<sup>9</sup>


Episode cost without reoperation (US)	Episode cost with reoperation (US)
<b>\$55,701</b>	<b>\$93,730</b>



# MAJOR CAUSES OF REOPERATION

A metanalysis of randomized evidence<sup>10</sup> published in 2015 analysing intertrochanteric fracture showed:

- 5.5% (248/4506) All-Cause Reoperation %**
- 3.5% Cut-Out %**
- 1.2% Distal Femur Fracture %**



# UNMET NEED



# SOLUTION

Cut-Out is the leading cause of failure of fixation in proximal femoral fractures accounting for half of all reoperations.

Cut-out rates for trochanteric hip nails with screw or helical blade femoral head fixation on average 3.5%.<sup>10,11</sup>



The X-Bolt is 25% stronger to push-out and 380% stronger in rotational stability vs hip screws or helical blades.<sup>12,13</sup>

The expanding wings give haptic feedback to the surgeon on the bone quality.

Large multi-center clinical trials with the X-Bolt have cut-out rates consistently <1%, before the recent design improvements giving better ergonomics and better tip-apex distance.<sup>14,15</sup>

Distal Fracture at the tip of either a short or long nail is the second leading cause of fixation failure at 1.2%.<sup>10,11</sup>

High stress concentration at the tip of a short or long nail, as well as mismatch of femoral curvature, lead to increased risk of fracture at distal nail tip.



The Pro-X1™ Nail has a distal taper and prongs to alleviate stress concentration at the tip of the nail. Biomechanical studies have shown this to significantly reduce strains in the femur at the nail tip.<sup>16</sup>

The Pro-X1™ long nails have been designed with a radius of curvature of 1.25m for nails 300mm to 375mm and 1.50m.

Over-engineered solutions contribute to excess operative time and deflect from the important clinical steps, such as fracture reduction, implant placement and tip-apex distance.<sup>11</sup>



The Metro Jig™ is a curved jig with flexi-drive that allows surgeons to 'operate around a corner' and facilitates faster surgery. It is especially useful in overweight and obese patients.

The simple surgical technique allows surgeons to concentrate on the important aspects such as fracture reduction, implant placement and tip-apex distance.

# FEMORAL HEAD CUT-OUT

## X-Bolt Design

- > Four orthogonal wings expand radially
- > Compacts surrounding cancellous bone
- > Tip-apex distance kept constant
- > No spinning of femoral head
- > Easily reversible

## Biomechanical Studies

Cyclic loading push-out testing showed the X-Bolt having a **28% greater peak force at cut-out**<sup>12</sup> versus a SHS screw, which has a thread outer diameter of 12.7mm. (Note: The thread outer diameters for Gamma3 and TFNA lag screws are smaller at 10.5mm and 10.35mm, respectively).




The femoral head rotates throughout a walking cycle or when standing from a seated position. Torque testing showed a **380% greater resistance to peak torque**<sup>13</sup> versus SHS screw and 87% greater than a helical blade.

X-Bolt placement is more forgiving than a screw. A study published in 2019 shows that push-out strength is superior for X-Bolt throughout a variety of tip-apex distances.<sup>17</sup>

## Clinical Studies

The X-Bolt has been the subject of two large randomised clinical trials in the UK, WHITE1<sup>14</sup> (n=100) and WHITE4<sup>15</sup> (n=1,128), where the X-Bolt XHS was compared to a generic SHS screw. Both constructs used a side plate and all other parameters were similar for both groups.

The overall **cut-out rate for the X-Bolt was 0.8%** (4/526). To date, there has not been a femoral head fixation device that has recorded a lower cut-out rate than the X-Bolt, and the clinical studies validate the superior biomechanical performance against cut-out.



X-Bolt Cut-Out Rate <sup>14,15</sup>	Traditional Screw Cut-Out Rate <sup>10</sup>
0.8%	3.5%





# DISTAL STRESS MODULATION

## Nail Design

- > Distal taper and prongs
- > Alleviates stress concentration at tip of nail

## Biomechanical Studies

Biomechanical studies have shown the design to consistently reduce strains in the femur at the nail tip, which thus will lessen the risk of a distal femoral fracture.<sup>16</sup>




X-BOLT NAIL      GAMMA3 NAIL

Change in strain vs intact %

# FASTER OPERATIVE TIME

**Average operating times** from the Norwegian Hip Fracture database<sup>5</sup> for each type of hip fracture operation are shown opposite. The mean operating time for a traditional short hip nail is 53 minutes. The average cost per minute for operating room use is estimated at \$36/minute.<sup>18</sup>

Just a 5-minute saving from use of the curved Metro Jig™ and the easy-to-use Pro-X1™ surgical technique will save a hospital \$180 per case and **\$36,000 annually** for a procedural volume of 200 cases per annum.




Mean Operating Times (minutes)



# ECONOMIC VALUE: THE X-BOLT® PRO-X1™ SOLUTION

## Reducing The High Cost of Reoperation

Reduction in reoperation from cut-out and distal femoral fracture will substantially reduce costs to the hospital and healthcare system. A sample budget impact analysis was developed to show the potential impact to a typical hospital that manages 200 hip fracture fixation per annum. The analysis evaluated the use of the X-Bolt Pro-X1™ nailing system and compared to a traditional screw-based hip nailing system, using data from published studies.<sup>10,14,15</sup> Results demonstrated that a typical hospital may recognize savings of up to **\$240,000** from the reduction in reoperations for fixation failure, mostly through a virtual elimination of distal femoral fracture and of 'cut-out'.



## Faster Hospital Discharge

Published clinical studies<sup>14,15</sup> have also shown better patient outcome scores for hip fracture patients treated with an X-Bolt device versus a screw-based device. Average difference in EQ-5D outcome score was +0.04 in favour of the X-Bolt at four months post-operatively. The difference seen was even greater at the 4-week mark at +0.09. Benchmarking against a similar hip fracture hemiarthroplasty trial<sup>19</sup> where a modern hemiarthroplasty achieved similar EQ-5D outcome score of +0.04, and a reduction a mean length of hospital stay of 0.67days (9.67days vs 9.00 days,  $p < 0.05$ ).

An estimated faster average hospital discharge of just 0.5 days translates into savings of \$1,186 per patient episode<sup>20</sup> in the US and **\$237,000** for a typical hospital with an annual procedural volume of 200 cases per year.



## Faster Surgical Time

The simplicity of the Pro-X1™ implant's design and use of the novel Metro Jig™ can save valuable minutes of the standard operating time. In particular the Metro Jig™ allows all-outside operation of the set screw deployment step and the jig decoupling step. The Metro Jig™ in effect allows the surgeon to 'operate around a corner' and is particularly useful in overweight or obese patients.

A 5-minute saving in average operating time results in savings of \$180 per case and **\$36,000** for a typical hospital<sup>18</sup> with an annual procedural volume of 200 cases per year.

## Total Hospital Savings

**up to \$513,000  
per annum**

## Total Savings

Combining the above, a typical hospital with an annual procedural volume of 200 cases per year could expect savings of up to **\$513,000** per annual using the X-Bolt® Pro-X1™.

# REFERENCES

1. Panula J, Pihlajamäki H, Mattila VM, Jaatinen P, Vahlberg T, Aarnio P, Kivelä SL. Mortality and cause of death in hip fracture patients aged 65 or older—a population-based study. *BMC musculoskeletal disorders*. 2011 Dec;12(1):1-6.
2. Blankart CR, van Gool K, Papanicolasl, Bernal-Delgado E, Bowden N, Estupiñán-Romero F, Gauld R, Knight H, AbionaO, Riley K, Schoenfeld AJ. International comparison of spending and utilization at the end of life for hip fracture patients. *Health Services Research*. 2021 Dec;56:1370-82.
3. Royal College of Physicians. National hip fracture database annual report 2017. <http://www.nhfd.co.uk/files/2017ReportFiles/NHFD-AnnualReport2017.pdf>
4. Authen AL, Dybvik E, FurnesO, GjertsenJE. Surgeon's experience level and risk of reoperation after hip fracture surgery: an observational study on 30,945 patients in the Norwegian Hip Fracture Register 2011–2015. *Acta orthopaedica*. 2018 Sep 3;89(5):496-502.
5. Norwegian National Advisory Unit on Arthroplasty and Hip Fractures. Helse Bergen HF, Department of Orthopaedic Surgery, HaukelandUniversity Hospital. <http://nrlweb.ihelse.net>. ISBN: 978-82-91847-25-2
6. American Society for Bone and Mineral Research. The Crisis. ASMBR: Secondary Fractures <https://www.secondaryfractures.org/the-crisis> (ASBMR, 2020).
7. Hernlund E, SvedbomA, IvergårdM, Compston J, Cooper C, StenmarkJ, McCloskey EV, Jönsson BK, KanisJA. Osteoporosis in the European Union: medical management, epidemiology and economic burden. *Archives of osteoporosis*. 2013 Dec;8(1):1-15.
8. Adeyemi A, DelhougneG. Incidence and economic burden of intertrochanteric fracture: a Medicare claims database analysis. *JBJS Open Access*. 2019 Mar 27;4(1).
9. Chitnis AS, VanderkarrM, Ruppenkamp J, Lerner J, Holy CE, Sparks C. Reoperations in intramedullary fixation of pertrochanteric hip fractures. *Journal of medical economics*. 2019 Jul 3;22(7):706-12.
10. Yu J, Zhang C, Li L, Kwong JS, XueL, Zeng X, Tang L, Li Y, Sun X. Internal fixation treatments for intertrochanteric fracture: a systematic review and meta-analysis of randomized evidence. *Scientific reports*. 2015 Dec 11;5(1):1-1.
11. HaidukewychGJ. Intertrochanteric fractures: ten tips to improve results. *JBJS*. 2009 Mar 1;91(3):712-9.
12. Gibson D, Keogh C, Morris S. A biomechanical study comparing the dynamic hip screw with an X-Bolt in an unstable intertrochanteric fracture model of the proximal femur. In *Orthopaedic Proceedings 2012 Sep (Vol. 94, No. SUPP\_XXXIX, pp. 164-164)*. The British Editorial Society of Bone & Joint Surgery.
13. Gosiewski JD, Holsgrove TP, Gill HS. The efficacy of rotational control designs in promoting torsional stability of hip fracture fixation. *Bone & joint research*. 2017 May;6(5):270-6.
14. Griffin XL, Parsons N, McArthur J, Achten J, Costa ML. The Warwick Hip Trauma Evaluation One: a randomised pilot trial comparing the X-Bolt Dynamic Hip Plating System with sliding hip screw fixation in complex extracapsular hip fractures: WHiTE (One). *The bone & joint journal*. 2016 May;98(5):686-9.
15. Griffin XL, Achten J, O'Connor HM, Cook JA, Costa ML, WHiTE Four Investigators. Effect on health-related quality of life of the X-Bolt dynamic plating system versus the sliding hip screw for the fixation of trochanteric fractures of the hip in adults: the WHiTE Four randomized clinical trial. *The Bone & Joint Journal*. 2021 Jan 4;103(2):256-63.
16. MacLeod AR, Rose H, Gill HS. Relative Fracture Risk for Plating and Nailing in Displaced Trochanteric Fractures. In *Orthopaedic Proceedings 2017 May (Vol. 99, No. SUPP\_9, pp. 72-72)*. The British Editorial Society of Bone & Joint Surgery.
17. Kahane S, Vaghela KR, Stammers J, Goldberg A, SmithamP. Biomechanical Study Comparing Cut-out Resistance of the XBolt® and Dynamic Hip Screw at Various Tip-Apex Distances. *Surgical technology international*. 2019 Nov 1;35:395-401.
18. Childers CP, Maggard-Gibbons M. Understanding costs of care in the operating room. *JAMA surgery*. 2018 Apr 1;153(4):e176233-.
19. Sims AL, Parsons N, Achten J, Griffin XL, Costa ML, Reed MR, CORNET Trainee Collaborative. A randomized controlled trial comparing the Thompson hemiarthroplasty with the Exeter polished tapered stem and Unitraxmodular head in the treatment of displaced intracapsular fractures of the hip: the WHiTE3: HEMI Trial. *Bone Joint J*. 2018 Mar;100(3):352-60.
20. Cost of inpatient day at hospitals in the U.S. 2019 by type. Frédéric Michas. Jan 17, 2022. <https://www.statista.com/statistics/630443/inpatient-day-hospital-costs-in-us-by-nonprofit-or-profit>

The information presented in this brochure is intended as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific X-Bolt products. Always refer to the package insert, product label and instructions for use before using any X-Bolt product. Surgeons must always rely on their own clinical judgement, training and expertise when deciding which products and techniques to use with their patients.

X-Bolt® is a registered trademark of X-Bolt Orthopaedics

X-Bolt Orthopaedics is a trading name of SOTA Orthopaedics Limited  
Registered in Ireland, No. 439651

European Patents: EP 2175790, EP 3496637, EP 2175790  
US Patents: US 9724141B2, US 8911446B2, US 11259854B2

VALUE-XBOLT-001.1 Date: 14 February 2022

For further information,  
or to order, please contact:

**X-Bolt Orthopaedics**  
Unit 5, Northwood Court  
Santry, Dublin D09AX54  
Ireland

**IE:** +353 1 443 3880

**US:** +1 813 358 3980

**W:** [www.x-bolt.com](http://www.x-bolt.com)